


DRY POWDER INHALER

DISKUS (6+ YEARS)



1
CLOSED



2
PUSH OPEN



3
SLIDE & CLICK



4
BREATHE OUT



5
DEEP BREATH IN & HOLD



6
RINSE & SPIT

IMPORTANT INFORMATION

Empty?: The window has numbers that count down the dose – “0” means it is empty. Throw it away.

Expired?: Check the expiry date on the back of the device.

Reminders:

- Keep the Diskus dry.
- Do not breathe into the device.
- Close after use.

Developed by the Community Pediatric Asthma Service, Calgary Health Region.

For more information on asthma medications and device demos, visit our website at:
www.calgaryhealthregion.ca/ican