

DRY POWDER INHALER


HANDIHALER



1 OPEN




2 LIFT



3 INSERT CAPSULE



4 CLOSE



5 PRESS & LET GO




6 BREATHE OUT



7 DEEP BREATH IN & HOLD



8 REMINDER



9 THROW AWAY

IMPORTANT INFORMATION

Expired?: Check the expiry date on the box or on the pack that holds the capsules.

- Reminder:**
- Keep the HandiHaler dry.
 - Do not breathe into the device.
 - Capsules are sensitive to light, moisture and extreme temperatures.
 - Make sure to remove only **ONE** capsule from the pack at a time.
 - If the dry powder gets on your hands, wash your hands well because the powder may irritate your eyes.

Developed by the Community Pediatric Asthma Service, Calgary Health Region.

For more information on asthma medications and device demos, visit our website at:

www.calgaryhealthregion.ca/ican