


DRY POWDER INHALER

DISKUS (6+ YEARS)



CLOSED




PUSH OPEN



SLIDE & CLICK



BREATHE OUT



DEEP BREATH IN & HOLD



RINSE & SPIT

IMPORTANT INFORMATION

Empty?: The window has numbers that count down the dose – “0” means it is empty. Throw it away.

Expired?: Check the expiry date on the back of the device.

- Reminders:**
- Keep the Diskus dry.
 - Do not breathe into the device.
 - Close after use.

Developed by the Community Pediatric Asthma Service, Calgary Health Region.

For more information on asthma medications and device demos, visit our website at:
www.calgaryhealthregion.ca/ican