



SOFT MIST INHALER

RESPIMAT - DAILY USE



1 KEEP CAP CLOSED
TURN BASE



2 FLIP OPEN CAP



3 BREATHE OUT



4 START SLOW BREATH,
PRESS BUTTON,
KEEP BREATHING IN



5 CLOSE CAP



6 REMINDER

IMPORTANT INFORMATION

Note: Place mouthpiece between your teeth and close your lips.

Empty?: The base has a dose tracker – “0” means it is empty. Base will lock. Throw it away.

Expired?: Check the expiry date on the device.

Reminders: If RespiMat isn't used for 1 week, **PRIME 1 TIME** (See Reverse Side).
If RespiMat isn't used for 3 weeks, **PRIME 3 TIMES** (See Reverse Side).
Do not get mist in eyes.

Cleaning: Wipe mouthpiece with a dry tissue or cloth.

This material is designed for information purposes only. It should not be used in place of medical advice, instruction and/or treatment. If you have specific questions, please consult your doctor or appropriate healthcare professional.

Illustrations developed by the Community Pediatric Asthma Service, www.ucalgary.ca/icancontrolasthma


For more information on inhaled medications and device demos, visit our website at:
www.ucalgary.ca/asthma




SOFT MIST INHALER

RESPIMAT - LOADING AND PRIMING


LOADING: REQUIRED ONCE, PRIOR TO FIRST USE



1
PRESS SAFETY
CATCH AND PULL
OFF CLEAR BASE.




2
INSERT NARROW
END OF CARTRIDGE
UNTIL IT CLICKS.
FIRMLY PRESS
AGAINST A HARD
SURFACE TO ENSURE
IT IS SECURELY
IN PLACE.



3
REPLACE
CLEAR BASE.

PRIMING CARTRIDGE: REQUIRED ON FIRST USE, FOLLOW T.O.P.




4
T
TURN
TURN BASE UNTIL
YOU HEAR "CLICK"



5
O
OPEN
FLIP OPEN CAP



6
P
PRESS
POINT DOWNWARD
AND PRESS BUTTON



7
R
REMINDER

- Close cap
- Repeat steps 4-6 until a cloud is visible
- Then repeat 3 more times before use

Note:

If Respimat isn't used for 1 week, **PRIME 1 TIME** (Steps 4-6).

If Respimat isn't used for 3 weeks, **PRIME 3 TIMES** (Steps 4-7).

Always keep cap closed when turning to avoid accidental activation of device.

Do not get mist in your eyes.

